

















File	Orifice Shaper 15/07	BlackJack 15/05	EZP 15/03	EZP 15/05	Shaping Files 03 Taper	Shaping Files 05 Taper	RT 20/06	RT 25/08
Length	17mm	17mm	21, 25, 31mm	21, 25, 31mm	21, 25, 31mm	21, 25, 31mm	21, 25mm	21, 25mm
Tip	Smart	Dynamic	Smart	Dynamic	Smart	Smart	Smart	Smart
Flutes	4	4	4	4	2	2	3	3
Handle	Silver	Black	Silver	Silver	Silver	Silver	Gold	Gold
Tip Color Band	ISO/White	ISO/White	ISO/White	ISO/White	ISO/20-45	ISO/20-45	ISO/Yellow	ISO/Red
Taper Color Band	Green	Red	Yellow	Red	Yellow	Red	Blue	Black
Stopper	Green	Red	Yellow	Red	Yellow	Red	Blue	Black
Cross Section								
Tip								
MFD	1.0mm	0.80mm	0.80mm	0.80mm	0.80mm	1.0mm	1.0mm	1.0mm
Speed	700-900	700-900	700-900	700-900	700-900	700-900	900-1000	900-1000
Torque	2.0	2.0	2.0	2.0	2.0	2.0	3.0	3.0

21mm Depth Marks	18, 19, 20
25mm Depth Marks	18, 19, 20, 22, 24
31mm Depth Marks	24, 25, 26, 28, 30

16mm Cutting Lgth for all but 15/05 + 15/07
12mm Cutting Lgth for 15/05 + 15/07

Maximum Flute Diameter

.03 Taper: 0.80mm
 .05 Taper (20-45): 1.0mm
 .05 Taper (15): 0.80mm
 - Including Blackjack -

.06 Taper: 1.0mm
 .07 Taper: 1.0mm
 .08 Taper: 1.0mm